LUNGClear™

CLARITY IN ONE BLOOD DRAW

Product Guidebook

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Product:
Lung Cancer Early Detection Blood Test

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Lung Cancer Early Detection

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2. Intended Use

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LUNGClear™ is a multi-marker laboratory developed diagnostic test that incorporates multiplex reverse transcriptase polymerase chain reaction (RT-PCR) and electrochemiluminescence immunoassay (ECLIA) intended for semi-quantitative detection of 12 lung neoplasia associated miRNAs and CEA protein biomarker from human platelet-poor plasma collected in BD Vacutainer® PPT™ Tubes (BD, NJ, USA). The LUNGClear test is performed on Applied Biosystems Veriti 96-well Thermocycler Instrument, QuantStudio Real-Time PCR Instrument with QuantStudio Software (Both instruments by Thermo Fisher Scientific Inc, MA, USA) and Cobas Pro Immunochemistry analyzer (Roche Diagnostics International Ltd, Rotkreuz, Switzerland) in M Diagnostics clinical laboratory.

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LUNGClear test is indicated for adults who are 40 years or older who are at average risk of lung cancer. The test is intended as an adjunctive tool complementary to other diagnostic imaging tests to aid in informing next steps in diagnosis of lung cancer.

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The results of LUNGClear test are presented as risk scores that must be interpreted by healthcare professionals alongside other laboratory findings as well as all relevant clinical evidence such as the patient's medical history. other laboratory findings. High-risk score may indicate the presence of lung cancer and should be followed up by diagnostic imaging tests and/or bronchoscopy. Low risk scores indicate low risk of developing lung cancer. Low risk scores should not be used as the sole basis for patient diagnoses and does not supplant clinical review by the physician. Individuals with a low-risk score is recommended to repeat blood test after 1 year or at an interval recommended by their physician.

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The LUNGClear test is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and in vitro diagnostic

procedures in a laboratory setting and under the supervision of a trained laboratory professional. The LUNGClear test is not intended for self-testing.

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3. Method of Detection

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1) Overview

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LUNGClear utilizes a proprietary RT-qPCR approach to detect multiple miRNA biomarkers associated with lung cancer. The test involves three steps:

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
- Step 1: RNA isolation from platelet-poor plasma samples
- Step 2: cDNA synthesis
- Step 3: Detection of miRNAs by quantitative PCR (qPCR)

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
The clinical lab workflow will start from the isolation of miRNA from platelet-poor plasma samples to cDNA synthesis, to detection of miRNAs by quantitative PCR (qPCR) before data analysis (amplification curve inspection).

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A risk score will then be generated incorporating CEA expression levels from platelet-poor plasma samples determined by ECLIA using our lab’s clinically validated algorithm.

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2) References

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- (1) Ying, L., Du, L., Zou, R., Shi, L., Zhang, N., Jin, J., ... & Su, D. (2020). Development of a serum miRNA panel for detection of early stage non-small cell lung cancer. *Proceedings of the National Academy of Sciences*, 117(40), 25036-25042.
- (2) Cheong, J. K., Tang, Y. C., Zhou, L., Cheng, H., & Too, H. P. (2022). Advances in quantifying circulatory microRNA for early disease detection. *Current opinion in biotechnology*, 74, 256-262.
- (3) Wan, G., Lim, Q. E., & Too, H. P. (2010). High-performance quantification of mature microRNAs by real-time RT-PCR using deoxyuridine-incorporated oligonucleotides and hemi-nested primers. *Rna*, 16(7), 1436-1445.
- (4) Ho, Y. K., Xu, W. T., & Too, H. P. (2013). Direct quantification of mRNA and miRNA from cell lysates using reverse transcription real time PCR: a multidimensional analysis of the performance of reagents and workflows. *PLoS one*, 8(9), e72463.

4. **Risk Score interpretation**

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Using RT-qPCR to detect 12 miRNA biomarkers associated with lung cancer and ECLIA to determine CEA expression levels, the risk score is calculated using a clinically validated algorithm. A high score obtained from the test may indicate the presence of lung cancer and hence the patient should undergo further evaluation with other diagnostic modalities such as LDCT, bronchoscopy and biopsy as required, in accordance with recognized clinical guidelines.

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5. High priority and urgent medical report situations

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In the event of a high-risk score, immediately contact the healthcare provider to inform them on the recommended follow up for the patient.

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6. Test Turn Around Time (TAT)

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7 - 10 Working days from sample receipt in lab to issuance of test report

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7. Sample collection requirements, equipment required, sample volume requirements

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1) Sample collection requirements

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- Using aseptic technique, obtain venous blood by using needle of 21G or 22G.

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NOTE: Do not apply the tourniquet for longer than 1 min as this may lead to hemolysis. Needles equal or finer than 23G are not recommended as it may lead to hemolysis as well.

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- Draw whole blood into BD Vacutainer® Serum Tube. The Serum tube is to be discarded; only around 1 ml of blood needs to be collected.

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- Gently invert the PPT tube 5 times to mix the blood with the tube additives. Be sure that tubes are NOT being shaken vigorously, as this can lead to a hemolyzed sample.

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- For plasma separation, centrifuge the blood tube at 1500 rcf for 15 min at room temperature using horizontal rotor (swing-out head).

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NOTE: Excessive centrifuge speed (over 2000 RCF) may cause tube breakage and exposure to blood and possible injury.

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6. After spinning, carefully aspirate 2.0 ml of the supernatant by using a 1 ml pipette to transfer 2 times (1.0 ml x 2 aspiration) and dispense into 15 ml centrifuge tubes.

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NOTE: Be very careful not to pick up red blood cells. Make sure to always aspirate from the top of the supernatant and leave at least 0.5 cm layer above the blood cells. Reduce aspiration volume if the volume of the supernatant is insufficient.

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7. Re-centrifuge the collected plasma in 15 ml centrifuge tubes at 2500 RCF at room temperature for 15 min.

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8. Transfer at least 1.5 ml of the supernatant into cryogenic vials and label as platelet-poor-plasma (PPP). Make sure to aspirate from the top of the supernatant and do not disturb any cells at the bottom of the tube. Leave at least 1/3 of total volume from the bottom of the tube. The bottom of the tube contains platelets which will confound the outcome of the Assay.

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9. Aliquot the PPP in appropriate volume (recommended 500 µl). Place all aliquots upright in a specimen box or rack in a -80°C or below freezer. All specimens should remain at -80°C or colder prior to shipping. The samples should not be thawed prior to shipping.

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10. All samples should be assessed for hemolysis by laboratory personnel with normal color vision. Refer to the hemolysis chart (Appendix 1). Re-collect blood if there is haemolysis. Store and transport the plasma at -80°C or dry ice until testing at the molecular diagnostics lab. Please fill up the Test Requisition Form (TRF) (Appendix 2).

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2) Equipment required

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- (1) BD Vacutainer needle 21G to 22G (Other manufacturer's 21 or 22G needle can be used.)

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- (2) BD Vacutainer one-use Holder (Cat. no.: 364815) (Other Manufacturer's holder can be used.)

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- (4) Centrifuge that can spin blood tubes

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Place barcode labels with two patient identifiers (Appendix 3) on the cryovials, as well as one barcode label on the TRF. If labels are unavailable, please use permanent marker to write down the same information directly on both cryovials and TRF.

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8. Sample Storage and Shipment

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- 1) Sample storage requirements

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-80°C for 6 months, -20°C for 3 days for processed plasma samples.

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- 2) Materials required

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- (1) Test Requisition Form (TRF)

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- (2) 3 cryovials per test with 2 patient identifiers

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- 3) Direct shipment protocol

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Before sending samples for shipment, please contact MiRXES to order courier for pick up.

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4) Storage/ Shipment protocol

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Please fold and insert the TRF into the Biohazard ziplock bag. Place cryovials in a cryobox and the box in the biohazard ziplock bag and keep frozen throughout. Ensure that the TRF and cryovials are adequately labelled with the relevant information. Note that the samples should not be thawed past this point.

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Please fill up the Test Requisition Form (TRF) (Appendix 3) with the following information:

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(6) Test information

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(7) Sample information

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(8) Blood taken date & time (DD/MM/YYYY __: __ 24-hour clock time)

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Please fill up the sample bar code label (Appendix 3) with the following information:

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(2) Medical Record Number

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11. Service lab contact (sample recipient)

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Please send your samples to our service lab M Diagnostics:

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Address: 10F, Hibiya Fort Tower, 1-1-1 Nishi-Shimbashi, Minato City, Tokyo, Japan

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電子メール : info.jp@mirxes.com

13. Report for test results

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Please note that the report for test results (Appendix 4) generated from M Diagnostics Lab in Singapore will be sent to clinics via digital email.

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Should there be any amendments to this "Product Guidebook" document, we will update and announce the change within a week's notice.

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15. Disclaimer

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Test performance may be altered based on factors such as excessive blood cell lysis before processing, sampling during treatment, and tissue heterogeneity. Possible sources of errors can result from trace contamination, rare technical errors, miRNA variants that interfere with analysis. This test is not intended to provide a definitive diagnosis of lung cancer and is not a substitute for imaging modalities (such as low-dose CT) and biopsy. Individuals with high-risk classification should be followed up by a healthcare professional.

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The interpretation and clinical annotations provided are intended solely for use by a medical professional and do not constitute medical advice by M Diagnostics. The clinical recommendation within the report is at the sole discretion of the physician. The clinical recommendation does not take into consideration of other risk factors such as genetic predispositions, other medical conditions, drug interactions, environmental and lifestyle factors. Decisions regarding prescriptions and treatment should not be solely based on this test and the information contained in this report.

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
The lung cancer risk classification may be updated over time with on-going clinical studies and trials. Information contained within this report is dependent on the date of generation and/or database version used to generate the report.

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Appendix 1. Haemolysis chart

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Circulating cell free miRNA are susceptible to contamination by blood cell miRNAs in the event of haemolysis.

Colorimetric assessment can provide an indication of the extend of haemolysis.

Re-collection of blood is recommended if the sample is hemolysed.

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Appendix 4. Sample Report

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 M Diagnostics <small>POWERED BY MIRXES</small>		30 Biopolis Street #09-05/06 Matrix Tel.: (65) 69503218 Company Registration No.: 202016378E	
<h2>LABORATORY TEST REPORT</h2> <h3>検査結果報告書</h3>			
Dr.	Name: 氏名	LUNG CLEAR DEMO	
	NRIC/Foreign ID: 国民/外国人登録番号	-	Sample Received Date: 09 Jan 2023 18:22 検体受領日
	Passport No./Others: パスポート等番号	985416	Order No.: M2301004908 発注番号
	Age/Gender: 年齢/性別	-/F	Accession No.: 220002737 検査受領番号
	DOB: 生年月日	Not Provided	Clinical Ref: AGE: 20
<h3>ONCOLOGY</h3>			
Collection Date: 検体採取日	09 Jan 2023		
Time: 検体採取時間	18:23:11		
TEST NAME 検査名	RESULT 検査結果	FLAG フラグ	UNIT 単位
LUNGClear (LC)	Low Risk		REFERENCE RANGE 基準範囲
LC Risk Score	12.7		
	Risk Classification リスク分類 Category カテゴリー <=12.7 - Low Risk <=12.7 - 低リスク Low risk of having lung cancer. The individual is recommended to repeat blood test after 1 year or at an interval recommended by a physician. 肺がんの低リスク。1年後、または医師が推奨する間隔で、血液検査を繰り返すことが推奨されます。		
	>12.7 - High Risk >12.7 - 高リスク High risk of having lung cancer. The individual is recommended to visit a specialist and consider a follow up examination. 肺がんの高リスク。専門医を受診し、さらなる検査を検討することが推奨されます。		
This report was validated and generated through the laboratory information system. No signature is required. Printed on 09 Jan 2023 18:24:32 by labadmin **Remarks: This report is only valid for clinical research/validation purposes.** Dr. Norman Walford, Medical Director, Ministry of Health, Singapore Laboratory License: L/210282/CLB/003/221			
			Page 1 of 3

mirxes

TO KNOW. TO ACT.

MiRXES Japan Co., Ltd.

10F, Hibiya Fort Tower, 1-1-1 Nishishimbashi, Minato City
Tokyo 105-0003 Japan

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